ROBOCAM
00-188 Warszawa
Dubois 10/5
Tel.: +48 22 894 59 92
www.robocam.info

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ROBOCAM
00-188 Warszawa
Dubois 10/5
Tel.: 22 894 59 92
www.robocam.info

Area of Application

Zirro disc are blanks from which items such as crowns and bridges in full shape, abutments, primary components or bar superstructures for dentures can be fabricated.

Safety first

Please read the manual carefully before you remove the zirconium oxide blank from the packaging. This advise if for your own protection when processing Zirro disc as well as the safety of your patients.

Warning! Processing zirconium disc and also Zirro disc produces dust which can lead to irritation of the eyes, skin and can damage the lungs. Always ensure that the extraction unit on your milling machine is in good working order. Wear goggles and a dust protection mask.

Zirro disc are produced and checked to the highest quality standards. It is absolutely essential to observe all the processing instructions in order to pass on this quality to patient. This will avoid any loss quality of the dental restoration in terms of aesthetics, fit or durability.

General handling instructions

Zirro disc are supplied in a party sintered state. At this stage they do not have the renowned, outstanding properties which characterise the final dental restoration. It is therefore necessary to handle them with care.

Please check the good immedietly after receipt for:

- damage to the packaging
- damage to the product (there must be no visible crack, material broken off or discoloration)

Store Zirro disc ideally:

- the original packaging
- in a dry place
- between 5 – 50 °C (40-120 F) temperature

When handling make sure that the Zirro disc:

- is not subjected to knocks or vibrations
- is not touched with wet hands
- does not come into contact with any sort of fluid (for example water, adhesives or inks)
- is not soiled


If any point in these instructions is not correct, the Zirro disc must not be used for the production of dental restorations.

Instruction for dentists

High quality care and subsequent long term patient satisfaction require high-grade materials and professional preparation and processing. You will obtain the best possible patient satisfaction if you follow the recommendations given below:

- prepare tooth as advised with chamfer type. It is premissible to prepare oh tooth under various angles
- depth of cut at the preparation margin at least 0.6 mm
- remove 0.6 to 1.0 mm of occlusal/incisal tooth material
- preparation angle: 4°-6°
- for laterals no more than two pontics may be positioned between two abutments teeth.

The restoration can be cemented with conventional phosphate or glass ionomer cement. As an adhesive the product RelyX Unicem from 3M ESPE is recommended

Instruction for dental laboratories

When designing crowns, Bridges and primary components, the following minimum dimensions should be adhered to:

- minimum wall thickness for anterior: 0.4 mm
- minimum wall thickness for laterals/ posterior: 0.6 mm
- min. wall thickness for abutment teeth: 0.6 mm
- cross-section of connectors anterior 9.0 mm²
- cross-section of connectors for laterals/posterior: 16 mm²
- restorations should be designed in such a way as to support the veneering ceramic at the cusp area, so that these can be veneered with an almost uniform layer of ceramic.

Notes on design are given in our "Guidelines for the design and processing of Zirro substructures".

Milling the disc

If the Zirro disc has arrived in an acceptable condition it is then ready to be processed in a Zirro CNC machine. It is important to position the block properly. Screws or any other elements holding the should be screwed on steadily in order to avoid a situation in which one side of the block is tense and the other is loose. With careful positioning the discs can be processed efficiently.

Removing the construction from the block

It is important to prevent the cut out construction from being subjected to local pressures. That is why one should:

- Cut the construction out of the block by sand blasting it with 50-100 µm aluminum oxide under the pressure of 2 atm
- Make sure that the connectors which join the block with the construction are cut out steadily so that none of its parts will be left.
- Processing should be conducted in a place that is tidy, clean and free of metallic fillings or any other substances which can cause discolorations during sinterization.

Coloring

The degree of Zirro's translucency is so high that it allows the material to be used in full-shape techniques. When coloring, it is recommended to proceed in compliance with one's general knowledge of coloring liquids and the guidelines specified in their manuals. The liquids should be applied on the crown one layer at a time according to the desired final effect. It is important for the successive layers – starting from the area surrounding the tooth's cervix up to its incisal margin and chewing surface – to blend with each other coloristically so that they are not visibly distinguishable from each other. Submerging the construction in water is not advisable. Before applying a layer of another color make sure that the brush is clean.

After permeating and sinterization, Zirro should not be subjected to machine processing as it can permanently discolor the surface. Finally, the construction should be coated with glazing.

Coloring construction

The traditional framework in ceramics should be colored by the soaking of the liquid coloring. We recommend to use acid-free liquid colorants. The time soaking in the liquid should be in accordance with the manufacturer's fluids. You should know that every coloring reduces the zirconia material transparency! It is recommended that the initial drying process before sintering in the Robocam oven.

Final sinterization of the construction

Final sinterization is a process in which the construction is made from Zirro material is subjected to precise thermal processing. Only then the substructure achieves the desired endurance which is necessary in order to guarantee durability and clinical success. Upon being thermally processed, the substructure shrinks to its final size as determined by the factor specified on the packaging. Sinterization requires specialized furnaces and near perfect cleanliness. Any contact of the Zirro material with metallic fillings or any other substances can result in permanent discoloration which would make achieving the desired final color an impossibility. It is recommended to regularly inspect the cleanliness of the furnace's chamber as well as the quality of its heating parts.

Sinterization pressure – 1450°C (2642 F)

Sintering furnace setup

During sinterisation of the Zirro substructures the sintering furnace should be setup as follows:

- from idle temp. To 900°C (1652 F) the temperature should rise by 10°C (50 F) per minute,
- substructure should be kept in 900°C (1652 F) for 30 min.
- from 900°C (1652 F) to 1450°C (2642 F) temperature should rise by 4°C (40 F) per minute,
- substructure should be kept in 1450°C (2642 F) for 2 hours

Cooling down:

- from 1450°C (2642 F) to 900°C (1652 F) the temperature should drop by 4°C (40 F) per minute,
- from 900°C (1652 F) to idle the temperature should drop freely.

Trimming and finishing the disc

Finally sinterized substructures and full-contour constructions should be processed only when it is absolutely necessary. Such processing would require the use of diamond or diamond-coated tools. It would also require the substructures and constructions to be cooled with water in order to avoid the materials' overheating which would likely change their structure and result in cracks.

Always abide by these guidelines:

- the thinner the wall of the prosthetic complement, the more caution is advised during machining
- during processing, the applied pressure should be kept to an absolute minimum-only diamond or diamond-coated tools should be used. The tools should be in perfect shape as ones that are worn out generate additional heat
- if possible, use tools which have a grain size smaller than 100 µm
- sharp edges should be avoided and smoothed out if possible.
- Do not use diamond-coated rubbers

Material properties and forms supply

The properties of hard sintered Zirro are set out in the following tables:

Shade white	51.05 +/- 2.02
Density	6.08 g/cm³
Open porosity	0 %
Twardość w Vickersach	1250 HV
Vickers hardness	4052.5 MPa
Bending strength (4 point)	1100 +/- 72.3 MPa
Yang's modulus	200.5 +/- 2.2 GPa
Fracture strength	5 +/- 0.16 MPaxm^{-1/2}
CTE (25-500 °C)	10.62 x 10⁻⁶
X-ray	W 100 % widoczna

Component	Content wt [%]
Y₂O₃	5.2
Al₂O₃	0.05
Fe₂O₃	<0.01
SiO₂	<0.02
ZrO₂	94

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